

AD-A125 449

COMPUTER ASSISTED DATA STORAGE AND RETRIEVAL IN
MUTAGENICITY TESTING II T. (U) LETTRMAN ARMY INST OF
RESEARCH PRESIDIO OF SAN FRANCISCO CA
W W JEDERBERG ET AL. NOV 82 LAIR-82-42TN

1/1

UNCLASSIFIED

F/G 6/3


NL

END

FILMED

221

DTIC



MICROCOPY RESOLUTION TEST CHART
NATIONAL BUREAU OF STANDARDS-1963-A

UNCLASSIFIED

SECURITY CLASSIFICATION OF THIS PAGE (When Data Entered)

REPORT DOCUMENTATION PAGE		READ INSTRUCTIONS BEFORE COMPLETING FORM
1. REPORT NUMBER LAIR Technical Note No. 82-42TN	2. GOVT ACCESSION NO. AD-A125449	3. RECIPIENT'S CATALOG NUMBER
4. TITLE (and Subtitle) Computer Assisted Data Storage and Retrieval In Mutagenicity Testing. II. The Drosophila Melanogaster Sex-Linked Recessive Lethal Asssay.		5. TYPE OF REPORT & PERIOD COVERED Final Dec 1979 to Dec 1980
7. AUTHOR(s) W. W. Jederberg, MS, CPT N. R. Powers, PhD, CPT		6. PERFORMING ORG. REPORT NUMBER
9. PERFORMING ORGANIZATION NAME AND ADDRESS Toxicology Group, Division of Research Support Letterman Army Institute of Research Presidio of San Francisco, CA 94129		8. CONTRACT OR GRANT NUMBER(s)
11. CONTROLLING OFFICE NAME AND ADDRESS Toxicology Group, Division of Research Support Letterman Army Institute of Research Presidio of San Francisco, CA 94129		10. PROGRAM ELEMENT, PROJECT, TASK 3M107-0109K UNIT NUMBERS WU 203
14. MONITORING AGENCY NAME & ADDRESS (if different from Controlling Office)		12. REPORT DATE November 1982
		13. NUMBER OF PAGES 26
		15. SECURITY CLASS. (of this report) UNCLASSIFIED
		15a. DECLASSIFICATION/DOWNGRADING SCHEDULE
16. DISTRIBUTION STATEMENT (of this Report) THIS DOCUMENT HAS BEEN APPROVED FOR PUBLIC RELEASE AND SALE, ITS DISTRIBUTION IS UNLIMITED.		
17. DISTRIBUTION STATEMENT (of the abstract entered in Block 20, if different from Report)		
18. SUPPLEMENTARY NOTES		
19. KEY WORDS (Continue on reverse side if necessary and identify by block number) Drosophila melanogaster Sex-Linked Recessive Lethal Assay Computer Assisted Data Storage and Retrieval		
20. ABSTRACT (Continue on reverse side if necessary and identify by block number) This report deals with the use of a computer for recording, storing and retrieving data generated from the Drosophila melanogaster Sex-Linked Recessive Lethal Assay.		

DTIC FILE COPY

DD FORM 1 JAN 73 1473

EDITION OF 1 NOV 65 IS OBSOLETE

UNCLASSIFIED

SECURITY CLASSIFICATION OF THIS PAGE (When Data Entered)

TECHNICAL NOTE NO. 82-42TN

COMPUTER ASSISTED DATA STORAGE AND RETRIEVAL
IN MUTAGENICITY TESTING

II. The *Drosophila melanogaster* Sex-linked Recessive Lethal Assay

WARREN W. JEDERBERG, MS, CPT MSC

and

NELSON R. POWERS, PhD, CPT MSC

TOXICOLOGY GROUP

DIVISION OF RESEARCH SUPPORT

NOVEMBER 1982

Toxicology Series 44

LETTERMAN ARMY INSTITUTE OF RESEARCH
PRESIDIO OF SAN FRANCISCO, CALIFORNIA 94129

88 AR 02 034

Computer-Assisted Data Storage and Retrieval in Mutagenicity Testing.
II. The Drosophila melanogaster sex-linked recessive lethal assay
--Jederberg and Powers


Reproduction of this document in whole or in part is prohibited except with the permission of the Commander, Letterman Army Institute of Research, Presidio of San Francisco, California 94129. However, the Defense Technical Information Center is authorized to reproduce the document for United States Government purposes.

Destroy this report when it is no longer needed. Do not return it to the originator.

Citation of trade names in this report does not constitute an official endorsement or approval of the use of such items.

In conducting the research described in this report, the investigation adhered to the "Guide for the Care and Use of Laboratory Animals," as promulgated by the Committee on Revision of the Guide for Laboratory Animal Facilities and Care, Institute of Laboratory Animal Resources, National Research Council.

This material has been reviewed by Letterman Army Institute of Research and there is no objection to its presentation and/or publication. The opinions or assertions contained herein are the private views of the author(s) and are not to be construed as official or as reflecting the views of the Department of the Army or the Department of Defense. (AR 360-5)


(Signature and date)

This document has been approved for public release and sale; its distribution is unlimited.

PREFACE

This technical note is the second in a series on the utilization of the computer facilities at Letterman Institute of Research to assist in mutagenicity testing as part of the institute's toxicology program. This report details the use of the computer in recording, storing and retrieving data. Its use in the Drosophila melanogaster sex-linked recessive lethal assay provides a form in which it is easy to view the raw data.

Accession For	
NTIS GRA&I	<input checked="checked" type="checkbox"/>
DTIC TAB	<input type="checkbox"/>
Unannounced	<input type="checkbox"/>
Justification	
By	
Distribution/	
Availability Codes	
Dist	Avail and/or Special
A	

DTIC
COPY
INSPECTED
2

ACKNOWLEDGMENTS

The authors wish to express their appreciation to the Information Sciences Group, Letterman Army Institute of Research, for making computer facilities available and to all those who have made suggestions for improving the utility of the program described herein.

COMPUTER-ASSISTED DATA STORAGE AND RETRIEVAL IN MUTAGENICITY TESTING.
II. The Drosophila melanogaster Sex-Linked Recessive Lethal Assay.

The mutagenicity testing of materials in compliance with federal regulations (4) requires several standards. Among these are the implementation of an extensive labelling system. A computer-assisted system has been established and described (2). Regulations also require storage of data. A computer program has been developed that records these new data from the Sex-Linked Recessive Lethal Assay (3), stores it in permanent computer files, and allows the user to see the data by means of a "print-out copy" as they are stored. The data also can be retrieved easily. Other programs which allow for a formatted print-out of a selected data set will also be described.

PROGRAM DESCRIPTION

The program DROSTOXDATA (Appendix A) records the raw data and stores it in the file DROSDATA. It records the number of failures, lethals and non-lethal offspring from each individual numbered male that was exposed to negative controls, positive controls, and test compounds. This recording is done for each of the four broods. In addition, the program makes a provision under the category of lethals so that the user may designate whether lethals are single or multiple. For each individual numbered male the program totals the number of failures, lethals (single or multiple), non-lethals, number of tests performed, and the mutation frequency. The program checks the total single and multiple lethals against the total as given by the individual broods. If they do not agree, an error message will appear and the user may re-enter the number of single and multiple lethals. If the total number of failures, lethals, and non-lethals for a given brood is greater than 99, only the brood data will be requested to be re-entered. As seen in the sample (Figure 1), the program displays the last run experiment number and the last male for which data were recorded.

The program requests the following information from the user (Figure 1):

(1) Select if more data from the same run are to be entered or if data for a new run are to be entered.

(2) Specify the run number if data for a new run are to be entered.

(3) Enter the identification number for each male fly tested. (The user is allowed up to 6 characters of any type.)

(4) Enter the identification of the compound tested. (The user is allowed up to 6 characters of any type.).

(5) Enter the data for the first brood as to the number of failures, lethals, and non-lethals. (If anything is greater than 99, an error message will be displayed and the user must re-enter the data for that brood. This procedure is repeated for broods 2 through 4).

If lethals are present:

(6) Enter the total number of single lethals.

(7) Then, enter the total number of multiple lethals. (The computer program will check the total number of single and multiple lethals against the total lethals in the broods. If there is an error, an error message will be displayed and the user will be requested to re-enter the number of single and multiple lethals. When total lethals agree, or there are none, the computer will type: TOTAL FAILURES, TOTAL LETHALS, TOTAL NON-LETHALS, TOTAL TESTS, and MUTATION RATE for data that have been entered.

(8) The user at this point may view how the data have been entered. If any of these data are incorrect the user may so indicate and return to Step 1.

(9) The user may then decide to continue by entering more data; beginning a new run or terminating the program.

(10) The user may log off the program and the data are stored in the DROSDATA file or the user may create a new file name to store the data.

(11) The user may receive a print-out of the data in their "raw" form (Figure 2), by giving the appropriate command to the operating system.

In addition, the user may receive a formatted print-out of the data accumulated by DROSTOXDATA by executing the program DROSRTPT (Figure 3). This formatted print-out (Figure 4) is generated from DROSRTPT and its associated subroutines RDDROSDATA, WRDROSDATA, and STRING (Appendices B, C, D and E).

The program, DROSRTPT, requests the following information from the

user (Figure 3):

(1) Type in the date; day, month, and year the program is being executed.

(2) Type in the name of the file to be printed out.

(3) Type in the appropriate command to the selected operating system name to receive a print-out of the data in the formatted form (Figure 4).

DISCUSSION

By utilizing these programs and the subroutines presented in this report, the raw data may be presented in a form which is easy to view and which saves time in analysis.

CONCLUSION

None.

RECOMMENDATION

None.

SAMPLE RUN OF X DROSTOXDATA

(Underlined entries are sample user inputs)

) X DROSTOXDATA

The computer will type:

LETTERMAN ARMY INSTITUTE OF RESEARCH
TOXICOLOGY GROUP: DROSOPHILA DATA

LAST RUN NUMBER RECORDED WAS: 37
LAST MALE DATA RECORDED WAS: T2-904

WOULD YOU LIKE:
TO ENTER MORE DATA FOR THE SAME RUN (= 1) ?
TO ENTER DATA FOR A NEW RUN (= 2) ?

ENTER CHOICE: (1 or 2)

1

MALE IDENTIFIER (AN/6) - ??

T2-905

TEST COMPOUND (AN/6) = ??

002MPT

ENTER FIRST BROOD DATA:

FAILURES (0-99), LETHALS (0-99), NONLETHALS (0-99):

1,0,24

ENTER SECOND BROOD DATA:

FAILURES (0-99), LETHALS (0-99), NONLETHALS (0-99):

0,0,25

ENTER THIRD BROOD DATA: FAILURES (0-99), LETHALS (0-99), NONLETHALS (0-99):

0,1,24

ENTER FOURTH BROOD DATA: FAILURES (0-99), LETHALS (0-99), NONLETHALS (0-99):

0,1,25

HOW MANY LETHALS WERE SINGLE ? (0-99)

1

Figure 1. Sample Run of X DROSTOXDATA

HOW MANY LETHALS WERE MULTIPLE ? (0-99)

0

*****ERROR*****

TOTAL LETHALS DON'T AGREE:

TOTAL LETHALS FROM BROODS = XX

TOTAL FROM SINGLES + MULTIPLES = XX

TOTAL FAILURES: 1 TOTAL LETHALS: 1

TOTAL NONLETHALS: 98 TOTAL TESTS: 99

MUTATION RATE FOR THESE DATA: 1.01

37 T2 905 002MPT 0,1,24 0,0,25 0,1,24 0,0,25 1 101 98 99 1.01

ARE THE DATA ALL RIGHT ? (1=YES, 0=NO)

1

WOULD YOU LIKE:

TO ENTER MORE DATA FOR THE SAME RUN (= 1)?

TO ENTER DATA FOR A NEW RUN (= 2) ?

TO ENTER NO MORE DATA (= 3) ?

ENTER CHOICE: 1, 2, OR 3. 1

)

Figure 1. Sample Run of X DROSTOXDATA

37 C1 908 1XFRUT 0 025 0 025 0 025 0 025 0 0 0 0 100 100 .00

37 C1 907 1XFRUT	0 025	0 025	0 025	0 025	0	0	0	0	100	100	.00
37 C1 908 1XFRUT	0 025	1 024	1 021	0 025	2	0	0	0	95	95	.00
37 C1 909 1XFRUT	1 024	0 025	0 124	0 025	1	1	0	1	98	99	1.01
37 C1 910 1XFRUT	1 024	1 024	0 025	0 025	2	0	0	0	98	98	.00
37 C1 911 1XFRUT	0 025	0 025	0 025	0 025	0	0	0	0	100	100	.00
37 C1 912 1XFRUT	0 025	0 025	1 024	0 025	1	0	0	0	99	99	.00
37 C1 913 1XFRUT	0 025	0 025	0 025	1 021	1	0	0	0	96	96	.00
37 P1 889 1EMST	2 419	1 222	0 520	0 0 0	3	11	011	61	72	15.28	
37 P1 890 1EMST	0 421	0 223	0 223	1 0 4	1	8	0 8	71	79	10.13	
37 P1 891 1EMST	0 223	0 124	1 222	0 0 0	1	5	0 5	69	74	6.76	
37 P1 892 1EMST	0 124	0 421	0 025	0 0 0	0	5	0 5	70	75	6.07	
37 P1 893 1EMST	0 421	0 421	0 124	0 111	0	10	010	77	87	11.49	
37 P1 894 1EMST	2 023	1 321	0 124	0 0 0	3	4	0 4	68	72	5.56	
37 P1 895 1EMST	1 321	3 022	8 215	0 0 0	12	5	0 5	58	63	7.94	
37 P1 896 1EMST	1 519	0 124	0 619	0 0 0	1	12	012	62	74	16.22	
37 P1 897 1EMST	0 223	2 122	1 321	0 0 0	3	6	0 6	66	72	8.33	
37 P1 898 1EMST	1 222	0 223	0 223	0 0 0	1	6	0 6	68	74	8.11	
37 T1 889 005MPT	0 025	0 025	0 025	0 025	0	0	0 0	100	100	.00	
37 T1 890 005MPT	0 025	0 025	1 024	0 025	1	0	0 0	99	99	.00	
37 T1 891 005MPT	0 025	0 025	0 025	0 025	0	0	0 0	100	100	.00	
37 T1 892 005MPT	0 025	0 025	0 025	0 025	0	0	0 0	100	100	.00	
37 T1 893 005MPT	0 025	0 025	1 024	1 0 0	1	0	0 0	99	99	.00	
37 T1 894 005MPT	0 025	3 121	0 025	0 025	3	1	0 1	96	97	1.03	
37 T1 895 005MPT	0 025	0 025	0 025	0 025	0	0	0 0	100	100	.00	
37 T1 896 005MPT	0 025	0 025	1 116	0 070	1	1	0 1	86	87	1.15	
37 T1 897 005MPT	0 025	0 025	1 024	0 016	1	0	0 0	90	90	.00	
37 T1 898 005MPT	1 024	0 025	0 025	0 025	1	0	0 0	99	99	.00	
37 T1 899 005MPT	0 025	0 025	0 025	0 020	0	0	0 0	95	95	.00	
37 T1 900 005MPT	1 024	0 025	0 025	0 025	1	0	0 0	99	99	.00	
37 T1 901 005MPT	0 025	0 025	0 025	0 025	0	0	0 0	100	100	.00	
37 T1 902 005MPT	0 025	2 023	0 025	0 025	2	0	0 0	98	98	.00	
37 T1 903 005MPT	1 024	0 0 0	0 0 0	0 0 0	1	0	0 0	24	24	.00	
37 T1 904 005MPT	0 025	0 025	0 025	1 024	1	0	0 0	99	99	.00	
37 T1 905 005MPT	1 024	0 025	1 024	0 025	2	0	0 0	98	98	.00	
37 T1 906 005MPT	0 025	0 025	0 025	0 025	0	0	0 0	100	100	.00	
37 T1 907 005MPT	0 025	0 025	0 025	1 024	1	0	0 0	99	99	.00	
37 T1 908 005MPT	0 025	1 024	0 025	1 016	2	0	0 0	92	92	.00	
37 T1 909 005MPT	0 025	0 0 0	0 0 0	0 0 0	0	0	0 0	25	25	.00	
37 T1 910 005MPT	0 025	0 025	0 025	0 025	0	0	0 0	100	100	.00	
37 T1 911 005MPT	1 024	0 025	0 025	0 025	1	0	0 0	99	99	.00	
37 T1 912 005MPT	2 023	0 025	0 025	0 025	2	0	0 0	98	98	.00	
37 T1 913 005MPT	3 022	0 025	0 0 0	0 0 0	3	0	0 0	47	47	.00	
37 T2 889 002MPT	5 020	0 025	0 025	0 025	5	0	0 0	95	95	.00	
37 T2 890 002MPT	0 025	0 025	0 025	0 0 9	0	0	0 0	84	84	.00	
37 T2 891 002MPT	0 0 0	0 0 0	0 0 0	0 0 0	0	0	0 0	0	0	.00	
37 T2 892 002MPT	0 015	0 025	0 025	1 024	1	0	0 0	89	89	.00	
37 T2 893 002MPT	0 020	0 025	0 025	0 025	0	0	0 0	95	95	.00	
37 T2 894 002MPT	1 024	14 011	0 025	0 025	15	0	0 0	85	85	.00	
37 T2 895 002MPT	0 025	1 024	0 010	0 0 0	1	0	0 0	59	59	.00	
37 T2 896 002MPT	0 025	0 0 0	0 0 0	0 0 0	0	0	0 0	25	25	.00	
37 T2 897 002MPT	1 024	0 0 0	0 0 0	0 025	1	0	0 0	49	49	.00	
37 T2 898 002MPT	1 024	0 025	0 025	0 025	1	0	0 0	99	99	.00	
37 T2 899 002MPT	0 025	0 025	0 025	0 025	0	0	0 0	100	100	.00	
37 T2 900 002MPT	0 025	0 025	0 025	0 0 0	0	0	0 0	13	13	.00	
37 T2 901 002MPT	0 025	0 025	0 025	1 024	1	0	0 0	94	94	.00	
37 T2 902 002MPT	0 025	0 025	0 025	0 025	0	0	0 0	100	100	.00	
37 T2 903 002MPT	0 025	0 025	0 025	0 025	3	0	0 0	97	97	.00	
37 T2 904 002MPT	0 025	0 025	0 025	0 025	0	0	0 0	100	100	.00	
37 T2 905 002MPT	1 024	0 025	0 024	0 025	1	1	0 1	98	99	1.11	
37 T2 906 002MPT	0 025	0 0 0	0 025	0 025	0	0	0 0	100	100	.00	
37 T2 907 002MPT	0 025	1 024	0 025	0 025	1	0	0 0	99	99	.00	

Figure 2. Raw Data of X DRUSTOXDATA

SAMPLE RUN OF X DROS RPT

(Underlined entries are sample user inputs)

) X DROS RPT

The computer will type:

ENTER DATA TO APPEAR ON REPORT (XXMONXX). 10Jun82

The computer will type:

NAME OF FILE TO BE REPORTED: DROSDATA

The computer will type:

REPORT READY IN 'SLRRPT'

)

Figure 3. Sample Run of X DROS RPT

LETTERMAN ARMY INSTITUTE OF RESEARCH
SEX-LINKED RECESSIVE LETHAL DROSOPHILA ASSAY
(RAA DATA PRINT-OUT)

PAGE: 30

DATE: 10JUN92

TEST	WAFER	COMPOUND	BLIND ALL FAIL NCL	BLIND B2 FAIL NCL	BLIND B1 FAIL NCL	BLIND B1 FAIL NCL	SUMMARY TOTALS TP NL ML TL TNL	TOTAL TESTS	MUTATION RATE (A)
37	12	896 002WPT	0 0 25	0 0 0	0 0 0	0 0 0	0 0 0 0 25	25	.00
37	12	897 002WPT	1 0 24	0 0 0	0 0 0	0 0 25	1 0 0 0 49	49	.00
37	12	898 002WPT	1 0 24	0 0 25	0 0 25	0 0 25	1 0 0 0 99	99	.00
37	12	899 002WPT	0 0 25	0 0 25	0 0 25	0 0 25	0 0 0 0 100	100	.00
37	12	900 002WPT	0 0 25	0 0 25	0 0 25	0 0 0 0	0 0 0 0 83	83	.00
37	12	901 002WPT	0 0 25	0 0 25	0 0 25	1 0 24	1 0 0 0 99	99	.00
37	12	902 002WPT	0 0 25	0 0 25	0 0 25	0 0 25	0 0 0 0 100	100	.00
37	12	903 002WPT	3 0 22	0 0 25	0 0 25	0 0 25	3 0 0 0 97	97	.00
37	12	904 002WPT	0 0 25	0 0 25	0 0 25	0 0 25	0 0 0 0 100	100	.00
37	12	905 002WPT	1 0 24	0 0 25	0 1 24	0 0 25	1 1 0 1 98	99	1.01
37	12	906 002WPT	0 0 25	0 0 25	0 0 25	0 0 25	0 0 0 0 100	100	.00
37	12	907 002WPT	0 0 25	1 0 24	0 0 25	0 0 25	1 0 0 0 99	99	.00
37	12	908 002WPT	1 0 24	0 0 25	0 0 15	0 0 7	1 0 0 0 71	71	.00
37	12	909 002WPT	0 0 25	0 0 25	0 0 25	0 0 25	0 0 0 0 100	100	.00
37	12	910 002WPT	0 0 25	0 0 25	0 0 25	0 0 17	0 0 0 0 92	92	.00
37	12	911 002WPT	0 0 25	1 0 24	0 0 25	0 0 9	1 0 0 0 83	83	.00
37	12	912 002WPT	0 0 25	0 0 25	0 0 25	0 0 20	0 0 0 0 95	95	.00
37	12	913 002WPT	0 0 25	0 1 24	0 0 25	1 0 21	1 1 0 1 95	96	1.04
38	01	914 102WPT	1 0 24	0 0 25	0 0 25	0 0 25	1 0 0 0 99	99	.00
38	01	915 102WPT	1 0 24	0 0 25	0 0 25	0 0 25	1 0 0 0 99	99	.00
38	01	916 102WPT	1 0 24	1 0 24	0 0 25	0 0 18	2 0 0 0 91	91	.00
38	01	917 102WPT	0 0 25	0 0 25	0 0 25	0 0 25	0 0 0 0 100	100	.00
38	01	918 102WPT	0 0 25	0 0 25	0 0 25	0 0 25	0 0 0 0 100	100	.00

EXPLANATORY NOTES: FATALITIES, LETHALS, NONLETHALS, MULTIPLE LETHALS, TL=TOTAL LETHALS,
TP=TOTAL FAILURES, SLENDING LETHALS, ML=MULTIPLE LETHALS, TL=TOTAL LETHALS,
TN=TOTAL NONLETHALS

Figure 4. Formatted Data of X DROSPT

REFERENCES

1. Food and Drug Administration. Good laboratory practices regulations. Federal Register 43(163):37336-37403, 1978.
2. JEDERBERG, W.W., R.A. WIRTZ and N.R. POWERS. Computer-Assisted Labelling in Mutagenicity Testing I. The Drosophila melanogaster Sex-Linked Recessive Lethal Assay. Technical Note No. 82-32TN, Letterman Institute of Research, Presidio of San Francisco, CA 94129, March 1982.
3. WIRTZ, R.A., N.R. POWERS and J.T. FRUIN. Mutagenicity Testing using the Drosophila melanogaster Sex-linked Recessive Lethal Assay. Institute Report No. 112, Letterman Army Institute of Research, Presidio of San Francisco, CA 94129, February 1982.

Appendix A.	Listing for DROSTOXDATA.FR	12
Appendix B.	Listing for DROSRPT.FR	17
Appendix C.	Listing for Subroutine RDDROSDATA.FR	19
Appendix D.	Listing for Subroutine WRDROSDATA.FR	20
Appendix E.	Listing for Subroutine STRING.FR	22

APPENDICES

Appendix A.

```

C --- PROGRAM NAME:  DROSTOXDATA.FR
C --- WRITTEN BY:    WARREN W JEDERBERG, CUTANEOUS HAZARDS
C --- INPUT:        AT TIME OF EXECUTION
C --- FORMAT FOR INPUT:  DASHER OR CRT
C --- SPECIAL EQUIPMENT: NONE
C --- PURPOSE:      TO STORE DATA FROM THE SEX-LINKED RECESSIVE LETHAL
C ---                DROSOPHILA TOXICOLOGICAL TEST, AND ALLOW CHECKING THE DATA
C ---                BEFORE IT IS STORED IN "DROSDATA" FILE
C ---
C --- REAL MR
C --- MR = MUTATION RATE FOR THE DATA
C --- INTEGER MNUM, SL, TL, TNL, TT, ML, CNAM, TTL
C --- DIMENSION CNAM(6),MNUM(6)
C --- MNUM = MALE INDENTIFICATION
C --- SL   = SINGLE LETHALS
C --- ML   = MULTIPLE LETHALS
C --- TL   = TOTAL LETHALS
C --- TNL  = TOTAL NONLETHALS
C --- TT   = TOTAL TESTS
C --- CNAM = COMPOUND INDENTIFIER
C ---
C --- INTEGER BFF,BFL,BFNL, BSF,BSL,BSNL, BTF,BTL,BTNL, BLF,BLL,BLNL
C --- BFF = FIRST BROOD FAILURES
C --- BFL = FIRST BROOD LETHALS
C --- BFNL = FIRST BROOD NONLETHALS
C ---
C --- BSF = SECOND BROOD FAILURES
C --- BSL = SECOND BROOD LETHALS
C --- BSNL = SECOND BROOD NONLETHALS
C ---
C --- BTF = THIRD BROOD FAILURES
C --- BTL = THIRD BROOD LETHALS
C --- BTNL = THIRD BROOD NONLETHALS
C ---
C --- BLF = FOURTH BROOD FAILURES
C --- BLL = FOURTH BROOD LETHALS
C --- BLNL = FOURTH BROOD NONLETHALS
C ---
C *** OUTPUT HEADER ***
99  WRITE (10,1)
1   FORMAT (///,10X,"LETTERMAN ARMY INSTITUTE OF RESEARCH",/,10X,
1      "TOXICOLOGY GROUP:  DROSOPHILA DATA",//)
C ---
C *** GET DATA LAST RECORDED ***
OPEN 2, "LDATA", ATT = "SIB"
READ (2,11) NR, MNUM
CLOSE 2

```

Appendix A.
(Continued)

```

11  FORMAT (I4,1X,6A1)
    WRITE (10,12) NR, MNUM
12  FORMAT (1X,"LAST RUN NUMBER RECORDED WAS:",2X,I4,/,
    11X,"LAST MALE DATA RECORDED WAS:",2X,6A1,/)
C ---
C *** OPTIONS (SAME RUN = 1, NEW RUN = 2) ***
301 CONTINUE
    WRITE (10,13)
13  FORMAT (10X," WOULD YOU LIKE:",/,15X,
    1"TO ENTER MORE DATA FOR THE SAME RUN (= 1) ?",/,15X,
    2"TO ENTER DATA FOR A NEW RUN (= 2) ?",/)
    ACCEPT " ENTER CHOICE: 1 OR 2 ",K
    IF (K.NE.1.AND.K.NE.2) GO TO 301
    IF (K.NE.2) GO TO 101
C ---
C *** GET NEW RUN NUMBER ***
999 CONTINUE
    TYPE
    ACCEPT " NEW RUN NUMBER ? ",NR
    IF (NR.LE.9999) GO TO 101
    TYPE
    TYPE " ++++++ERROR+++++"
    TYPE " RUN NUMBER > 9999"
    GO TO 999
101 CONTINUE
C ---
C *** GET MALE IDENTIFIER ***
    TYPE
    TYPE " MALE IDENTIFIER (AN/6)= ??"
    READ (11,8) (MNUM(K), K=1,6)
8   FORMAT (6A1)
C ---
C *** GET COMPOUND IDENTIFIER ***
    TYPE
    TYPE " TEST COMPOUND (AN/6)= ??"
    READ (11,7) (CNAM(K), K=1,6)
7   FORMAT (6A1)
C ---
C *** ENTER BROOD DATA
2   FORMAT (/, "ENTER FIRST BROOD DATA:")
3   FORMAT (/, "ENTER SECOND BROOD DATA:")
4   FORMAT (/, "ENTER THIRD BROOD DATA:")
5   FORMAT (/, "ENTER FOURTH BROOD DATA:")
6   FORMAT ("FAILURES (0-99), LETHALS (0-99), NONLETHALS (0-99):")
C ---

```

Appendix A.
(Continued)

```

C --- ENTER FIRST BROOD DATA
10  WRITE (10,2)
    WRITE (10,6)
    ACCEPT "**", BFF,BFL,BFNL
    IF (BFF.GT.99.OR.BFL.GT.99.OR.BFNL.GT.99) GO TO 10

C ---
C --- ENTER SECOND BROOD DATA
20  WRITE (10,3)
    WRITE (10,6)
    ACCEPT "**", BSF,BSL,BSNL
    IF (BSF.GT.99.OR.BSL.GT.99.OR.BSNL.GT.99) GO TO 20

C ---
C --- ENTER THIRD BROOD DATA
30  WRITE (10,4)
    WRITE (10,6)
    ACCEPT "**", BTF,BTL,BTNL
    IF (BTF.GT.99.OR.BTL.GT.99.OR.BTNL.GT.99) GO TO 30

C ---
C --- ENTER FOURTH BROOD DATA
40  WRITE (10,5)
    WRITE (10,6)
    ACCEPT "**", BLF,BLL,BLNL
    IF (BLF.GT.99.OR.BLL.GT.99.OR.BLNL.GT.99) GO TO 40

C ---
C *** CALCULATE TOTALS ***
    ML = 0.00
    SL = 0.00
    TL = BFL + BSL + BTL + BLL
    IF (TL.LE.0.1) GO TO 150
    TYPE
    TYPE
50  TYPE " HOW MANY LETHALS WERE SINGLE ? (0-99)"
    ACCEPT "**",SL
    TYPE " HOW MANY LETHALS WERE MULTIPLE ? (0-99)"
    ACCEPT "**",ML

C ---
C --- CHECK TOTAL LETHALS FROM BROODS AND FROM QUERY
    TTL = SL + ML
    IF (ABS(TTL-TL).LE.0.1) GO TO 150
    TYPE
    TYPE "*****ERROR*****"
    TYPE " TOTAL LETHALS DON'T AGREE:"
    TYPE " TOTAL LETHALS FROM BROODS = ",TL
    TYPE " TOTAL FROM SINGLES + MULTIPLES = ",TTL
    TYPE
    GO TO 50

```

Appendix A.
(Continued)

```

150 CONTINUE
    TNL = BFNL + BSNL + BTNL + BLNL
    TT = TNL + TL
    IF (TT.LE.100) GO TO 60
    TYPE "*****ERROR*****"
    TYPE "TOTAL TESTS ARE MORE THAN 100, RE-ENTER ALL DATA"
    GO TO 99
60 CONTINUE
    MR = TL*(100.00/TT)
    TF = BFF + BSF + BTF + BLF
C ---
C *** CHECK DATA ****
C --- DISPLAY CALCULATED DATA ***
    WRITE (10,70) TF,TL,TNL,TT,MR
70 FORMAT (/,"TOTAL FAILURES:",1X,I3,2X,"TOTAL LETHALS:",1X,I3,/,
1"TOTAL NONLETHALS:",1X,I3,2X,"TOTAL TESTS:",1X,I3,/,
2"MUTATION RATE FOR THESE DATA:",2X,F6.2)
C --- DATA AS WILL BE RECORDED
    WRITE (10,80)
80 FORMAT (//,"DATA WILL BE RECORDED AS:")
    WRITE (10,100) NR,MNUM,CNAM,BFF,BFL,BFNL,BSF,BSL,BSNL,
1BTF,BTL,BTNL,BLF,BLL,BLNL,TF,SL,ML,TL,TNL,TT,MR
    TYPE
81 TYPE "ARE THE DATA ALL RIGHT ? (1=YES,0=NO)"
    ACCEPT "**",I
    IF (I.NE.O.AND.I.NE.1) GO TO 81
    IF (I.NE.1) GO TO 99
C ---
C *** STORE DATA IN DROSDATA ***
    OPEN 1, "DROSDATA", ATT = "SA"
    WRITE (1,100) NR,MNUM,CNAM,BFF,BFL,BFNL,BSF,BSL,BSNL,
1BTF,BTL,BTNL,BLF,BLL,BLNL,TF,SL,ML,TL,TNL,TT,MR
100 FORMAT (I4,1X,6A1,1X,6A1,4(1X,3I2),1X,I2,1X,3I2,
12(1X,I3),1X,F6.2)
    CLOSE 1
C ---
C *** RECORD LAST RUN NUMBER AND MALE IN LDATA
    OPEN 2, "LDATA", ATT = "SO"
    WRITE (2,11) NR, MNUM
    CLOSE 2
C ---

```

Appendix A.
(Continued)

```
C *** OPTION FOR MORE DATA OR LOG OFF ***
200 CONTINUE
    WRITE (10,14)
14  FORMAT (//,10X,"WOULD YOU LIKE:",/,15X,
1    "TO ENTER MORE DATA FOR THE SAME RUN (= 1) ?",/,15X,
2    "TO ENTER DATA FOR A NEW RUN (= 2) ?",/,15X,
3    "TO ENTER NO MORE DATA (= 3) ?",/)
    ACCEPT " ENTER CHOICE: 1, 2, OR 3 ",K
    IF (K.NE.1.AND.K.NE.2.AND.K.NE.3) GO TO 200
    IF (K.NE.2.AND.K.NE.3) GO TO 101
    IF (K.NE.3) GO TO 999
C ---
    END
```

Appendix B.

```

C *** PROGRAM NAME:          DROSRPT.FR
C *** WRITTEN BY:           WARREN JEDERBERG
C *** PURPOSE:              TO GENERATE A REPORT OF THE
C ***                      RAW DATA AS STORED IN 'DROSDATA'

      INTEGER NR,MNUM,CNAM,BFF,BFL,BFNL,BSF,BSL,BSNL,
1BTF,BTL,BTNL,BLF,BLL,BLNL,TF,SL,ML,TL,TNL,TT,SFAM

      REAL MR

      COMMON/DREC/NR,MNUM(6),CNAM(6),BFF,BFL,BFNL,BSF,BSL,BSNL,
1BTF,BTL,BTNL,BLF,BLL,BLNL,TF,SL,ML,TL,TNL,TT,MR

      LOGICAL EOF

      DIMENSION NATE(5),SFAM(20)

C --- GET REPORT DATE

      TYPE "Enter Date to appear on report: (XX-MON-XX)"
      READ (11,400) (NATE(I),I=1,5)
400 FORMAT (5A2)

C --- INITIALIZE PAGE AND LINE
      PAGE = 0
      LINE = 0

C --- INITIALIZE EOF
      EOF = .FALSE.

C --- GET AND OPEN FILE TO BE REPORTED
      WRITE (10,300)
300 FORMAT (/"NAME OF FILE TO BE REPORTED: ",Z)
      READ (11,310) SFAM
310 FORMAT (20A2)
      CALL STRING (SFAM,20)
      OPEN 1, SFAM, ATT="SIB",ERR=100
      READ (1,1)
1  FORMAT (/)
      GO TO 120
100 CONTINUE
      WRITE (10,6) SFAM
6  FORMAT (//,"***** ERROR **** ERROR **** ERROR *****",/,
14X,"FILE ",20A2," NOT FOUND . . .")
      STOP

```

Appendix B.
(Continued)

```
C --- DELETE OLD REPORT
      DELETE "SLRRPT"

C --- OPEN REPORT FILE SLRRPT
      120 OPEN 2, "SLRRPT", ATT = "SOP"

C --- READ DATA RECORD
      140 CALL RDDROSDATA (EOF)
          IF(EOF) GO TO 220

C --- REPORT DATA

          CALL WRDROSDATA (LINE,PAGE,NATE(1))
          GO TO 140

C --- EOF DETECTED
      220 CLOSE 1
          IF (LINE.NE.O.OR.PAGE.NE.O) GO TO 250
          TYPE
          TYPE " EMPTY DATA FILE"
          CLOSE 2
          STOP

      250 CLOSE 2
          TYPE
          TYPE " REPORT READY IN 'SLRRPT'"
          END
```


Appendix C.

```
COMPILER NOSTACK
C *** PROGRAM NAME:          RDDROSDATA.FR
C *** WRITTEN BY:           WARREN JEDERBERG
C *** PURPOSE:              TO READ DATA FOR MAKING REPORTS
C ***                      FOR THE SLRL-DROSOPHILA ASSAY

SUBROUTINE RDDROSDATA (EOF)
COMMON/DREC/NR,MNUM(6),CNAM(6),BFF,BFL,BFNL,BSF,BSL,BSNL,
1BTF,BTL,BTNL,BLF,BLL,BLNL,TF,SL,ML,TL,TNL,TT,MR

INTEGER NR,MNUM,CNAM,BFF,BFL,BFNL,BSF,BSL,BSNL,BTF,BTL,BTNL,
1BLF,BLL,BLNL,TF,SL,ML,TL,TNL,TT

REAL MR

LOGICAL EOF

READ (1,100, END= 120) NR,MNUM,CNAM,BFF,BFL,BFNL,BSF,BSL,BSNL,
1BTF,BTL,BTNL,BLF,BLL,BLNL,TF,SL,ML,TL,TNL,TT,MR

100 FORMAT(I4,1X,6A1,1X,6A1,4(1X,3I2),1X,I2,1X,3I2,
12(1X,I3),1X,F6.2)
RETURN

120 EOF = .TRUE.
RETURN
END
```

Appendix D.

```
C *** PROGRAM NAME:          WRDROSDATA.FR
C *** WRITTEN BY:           WARREN JEDERBERG
C *** PURPOSE:             THIS SUBROUTINE HANDLES PAGING AND
C ***                     OUTPUTTING OF RECORDS FOR DROSAPT.
```

```
      SUBROUTINE WRDROSDATA (LINE,PAGE,IDATE)
      INTEGER NR,MNUM,CNAM,BFF,BFL,BFNL,BSF,BSL,BSNL,BTF,BTL,BTNL,
      1BLF,BLL,BLNL,TF,SL,ML,TL,TNL,TT
```

```
      DIMENSION IDATE(5)
      REAL MR
```

```
      COMMON/DREC/NR,MNUM(6),CNAM(6),BFF,BFL,BFNL,BSF,BSL,BSNL,BTF,BTL,
      1BTNL,BLF,BLL,BLNL,TF,SL,ML,TL,TNL,TT,MR
      IF (LINE.GT.0.AND.LINE.LT.54) GO TO 120
```

```
      PAGE=PAGE + 1
      WRITE (2,100) (IDATE(I),I=1,5), PAGE
```

```
100 FORMAT (1H1,42X,"LETTERMAN ARMY INSTITUTE OF RESEARCH",
      1/,39X,"SEX-LINKED RECESSIVE LETHAL DROSOPHILA ASSAY",
      2/,50X,"(RAW DATA PRINT-OUT)",
      3//," DATE:",2X,5A2,94X,"PAGE:",1X,I3,
      4//,27X,"BROOD #1:   BROOD #2:   BROOD #3:   BROOD #4:",9X,
      5"SUMMARY TOTALS",4X,"TOTAL",4X,"MUTATION",/,
      61H+,26X,4(9(" "),4X),3X,19(" "),
      7//," RUN:",2X,"MALE:",2X,"COMPOUND:",4X,4("F: L: NL:",4X),3X,
      82X,"TF SL ML TL TNL",3X,"TESTS:",3X,"RATE (%)")
```

```
      LINE = 3
```

```
120 WRITE (2,140) NR,MNUM,CNAM,BFF,BFL,BFNL,BSF,BSL,BSNL,
      1BTF,BTL,BTNL,BLF,BLL,BLNL,TF,SL,ML,TL,TNL,TT,MR
```

```
140 FORMAT (/ ,I4,3X,6A1,1X,6A1,7X,4(3(I2,1X),4X),5X,4(I2,1X),I3,
      14X,I3,6X,F6.2)
      LINE = LINE + 2
```

```
      IF (LINE.GE.54) GO TO 150
      RETURN
```

```
150 WRITE (2,160)
```

Appendix D.
(Continued)

```
160 FORMAT (/, " EXPLANATORY NOTES: F=FAILURES, L=LETHALS, NL=NONLETHALS", /,  
120X, "TF=TOTAL FAILURES, SL=SINGLE LETHALS, ML=MULTIPLE LETHALS, TL=",  
2"TOTAL LETHALS, ", /, 20X, "TNL=TOTAL NONLETHALS", /)
```

```
LINE=0  
RETURN
```

```
END
```

Appendix E.

```
C --- TO CONVERT FILE NAMES TO USE IN OPEN STATEMENTS
      SUBROUTINE STRING(LINE,LLEN)
      COMPILER STATIC
      DIMENSION LINE (LLEN)
      DO 100 I=1,LLEN
      IF (LINE(I).EQ." ") LINE(I)=0
      IF (FLD(LINE(I),9,16).EQ.FLD(" ",9,16))FLD(LINE(I),9,16)=0
100  CONTINUE
      RETURN
      END
```

END

FILMED

4-83

DTIC